M MILLIPORE

La ciencia del Western Blot

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D. en C. Samuel García Nieto
Asesor Life Sciences, Millipore

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Western Blot MILLIPORE

Gel Electrophoresis ⇒ **WESTERN BLOTTING** ⇒ **DETECTION**

Separates a complex mixture of proteins into individual components. Transfers separated proteins to a microporous membrane, e.g., PVDF. Identifies presence or absence of specific proteins

- Sample prep chemicals
- Electrophoresis hardware
- Pre-cast gels
- Protein loading buffers
- Running buffers
- Gel stains
- Protein standards

- Blotting (tank or semi-dry) hardware
- Transfer buffers
- NC or PVDF membranes
- Blotting paper
- Membrane sandwiches
- Blot stains

- Blocking reagents
- Wash buffers
- Primary antibody
- Secondary antibody
- Detection reagents
 - Chromogenic
 - Chemiluminescent
 - Fluorescent
- X-ray film/processing equipment
- Imaging Equipment

Preparar muestras MILLIPORE

Sample Prep → Gel Electro-phoresis → Membrane Transfer → Blocking → Antibody Addition → Detection

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Sample in Solution

Electroforesis & Transferencia MILLIPORE

Sample Prep → Gel Electro-phoresis → Membrane Transfer → Blocking → Antibody Addition → Detection

❖ Separar proteínas por electroforesis en gel

- 1-D
- 2-D

Cell Lysate or Protein mixture

Electrophoresis

SDS-Polyacrylamide Gel

Electroforesis & Transferencia MILLIPORE

Sample Prep → Gel Electrophoresis → Membrane Transfer → Blocking → Antibody Addition → Detection

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The diagram shows two systems for protein transfer. The **Wet Tank System** consists of a gel, a membrane, and a nitrocellulose membrane, with a buffer reservoir. The **Semi-Dry Tank System** consists of a gel, a membrane, and a nitrocellulose membrane, with a buffer reservoir and a support layer.

Membranas para blots MILLIPORE

- ❖ Polyvinylidene fluoride (PVDF)
 - 0.45 μm y 0.22 μm
 - Baja fluorescencia
- ❖ Nitrocelulosa
 - 0.45 μm
- ❖ Nylon (Neutral y Cargado positivamente)
 - 0.45 μm

The image shows a microscopic view of a membrane surface with a 3D protein model overlaid, illustrating the interaction between the protein and the membrane.

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Problema #1

Transferencia pobre

ADVANCING LIFE SCIENCE TOGETHER™
Research. Development. Production.

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MILLIPORE

Solución #1

- ❖ Ajustar condiciones de transferencia
 - Metanol
 - SDS
 - Voltaje
- ❖ Equilibrar gel y membrana en buffer de transferencia
 - Incubar por 5 - 30 mins
- ❖ Buffers:
 - Towbin buffer (Tris/glycine transfer buffer): 25 mM Tris base, 192 mM glycine, 10% (v/v) methanol, pH 8.3)
- ❖ Original Reference:

→ Towbin H, Staehelin J, Gordon J. "Electrophoretic transfer of proteins

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MILLIPORE

Solución #1


Figure 9. Effect of equilibration time on electrotransfer of BSA to Immobilon-P transfer membrane. ¹²⁵I-labeled BSA was resolved by SDS-PAGE on a 10 – 20% gradient gel. After equilibration for the times noted, the BSA was transferred to Immobilon-P transfer membrane, backed up with Immobilon-PSQ transfer membrane, in a tank transfer system using 25 mM Tris, 192 mM glycine, and 10% methanol, as the transfer buffer. The system was run at 8 V/cm interelectrode distance. At the end of the 2-hour transfer period, the gel and membranes were stained. The BSA bands were excised and counted.

M MILLIPORE

Problema #2

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ADVANCING LIFE SCIENCE TOGETHER™
Research. Development. Production.

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MILLIPORE

0.2um PVDF retiene más que 0.45um

Prolonged electrotransfer of proteins using Immobilon-P and Immobilon-P92 transfer membranes. Molecular weight standards [lanes 1,3,5,7] and calf liver lysate [lanes 2,4,6,8] were transferred to Immobilon-P or Immobilon-P92 membranes by the tank transfer method and stained with Coomassie Blue. A sheet of Immobilon-P92 transfer membrane was placed behind the primary membranes to capture proteins that passed through them. [Lanes 5 and 6 behind Immobilon-P; lanes 7 and 8 behind Immobilon-P92.]

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Sample Prep → Gel Electro-phoresis → Membrane Transfer → Blocking → Antibody Addition → Detection

Transfer Verification

A. Ponceau-S red B. CPTS C. Coomassie- Brilliant Blue

M MILLIPORE

Problema #3

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ADVANCING LIFE SCIENCE TOGETHER™
Research. Development. Production.

MILLIPORE

Western Blotting Workflow

Sample Prep → Gel Electro-phoresis → Membrane Transfer → Blocking → Antibody Addition → Detection

The diagram illustrates the Western Blotting Workflow. It starts with a gel containing protein bands. The process involves:

- Sample Prep**: Preparing the protein samples.
- Gel Electro-phoresis**: Running the gel to separate proteins by size.
- Membrane Transfer**: Transferring the protein bands from the gel to a membrane.
- Blocking**: Treating the membrane to prevent non-specific binding.
- Antibody Addition**: Adding a primary antibody (Add 1° Ab) that binds to the target protein, followed by a secondary antibody (Add 2° Ab) that binds to the primary antibody.
- Detection**: Visualizing the bound antibodies, often using a substrate that produces a color change or light signal.

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The diagram shows the structure of an IgM antibody. It is a pentamer, consisting of five Y-shaped units. Each unit has two heavy chains (labeled 'heavy chain') and two light chains (labeled 'light chain'). The variable regions are labeled as 'VH domain' and 'VL domain'. The constant region is labeled as 'C1'.

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The diagram compares the structures of polyclonal and monoclonal antibodies. Polyclonal antibodies are represented by a mixture of different Y-shaped molecules, each recognizing a different epitope. Monoclonal antibodies are represented by a single, uniform Y-shaped molecule, recognizing only one specific epitope.

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
Immobilon-P

Immobilon-PS^Q

Std Method	Extra TBST wash	High Salt wash

MILLIPORE

Problema #6

Dificultades al re-analizar

ADVANCING LIFE SCIENCE TOGETHER™
Research. Development. Production.

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MILLIPORE

Before ReBlot™ Plus	After ReBlot™ Plus	Reprobe

Before ReBlot™ Plus: Western blot probed with rabbit antibodies to peptide from angiotensin II type 1 receptor, 20 second exposure.
After ReBlot™ Plus: Western blot stripped with ReBlot™ Plus Mild stripping solution.
Reprobe: Western blot reprobed with rabbit antibodies to peptide from glutamate receptor delta 1/2, 20 second exposure.

M **MILLIPORE**

Problema #7

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ADVANCING LIFE SCIENCE TOGETHER™
Research. Development. Production.

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M **MILLIPORE**

AYUDA O DUDAS

http://www.millipore.com/immunodetection/id3/western_blotting

Western Blotting

Factors influencing immunodetection

- Buffers
- Blocking
- Antibodies
- Washing
- Double labeling
- Detector Substrates
- Reporting Immovable PVDF Transfer Membranes

Since its introduction in 1979 (Towbin et al., 1979), protein blotting has become a mainstay of research laboratories. Its technology used to identify primary or proteins in complex samples or to monitor protein expression and purification. The primary protein blotting procedure involves an SDS-PAGE or 2D-PAGE gel, which is then transferred onto a nitrocellulose membrane. This blotting method may provide valuable information about total protein expression levels and can be performed in either a single or multiple steps. For more information on protein blotting, visit www.millipore.com/immunodetection/id3/western_blotting.

Product News & Releases

New Product: Quanta Millipore antibodies and kits. They are intended for immunoblotting and protein research applications.

New Quanta™ Millipore antibodies.

M **MILLIPORE**

AYUDA O DUDAS

samuel_nieto@millipore.com

ADVANCING LIFE SCIENCE TOGETHER™
Research. Development. Production.